



Pinus radiata
MultitAAG[®]
GENOTYPING MADE EASY



One-tube reaction
Pinus radiata
genotyping kit & service

Pinus radiata
MULTITAAG[®]

→ **Get a fully informative fingerprint
in one reaction**

Fingerprint specifications

Dye	Locus	Repeat motif	Product size range (bp)	Exclusion power
6-FAM	TPRA01	(CT) _n	92-162	0,451
6-FAM	TPRA03	(AC) _n (AT) _n	243-330	0,735
6-FAM	TPRA04	(CA) _n	335-380	0,256
HEX	TPRA06	(CT) _n (CA) _n	193-239	0,715
HEX	TPRA07	(TA) _n (TG) _n (ATGT) _n	256-368	0,451
HEX	TPRA08	(AG) _n	382-430	0,476
HEX	TPRA09	(TG) _n	444-473	0,599
NED	TPRA10	(TC) _n	104-204	0,258
NED	TPRA11	(GA) _n	210-242	0,383
NED	TPRA12	(CA) _n (TA) _n	282-374	0,451
NED	TPRA13	(TG) _n	403-485	0,124

→ **Choose the most convenient option**

Option 1:
Purchase our
MultiTAAG[®] kit

Option 2:
Use our full
genotyping service

For more information,
please contact us at:
forestry@taag-genetics.com
www.taag-genetics.com